



Aktin Chemicals, Inc.

No.52, Yongfeng Road, Hi-tech Zone, Chengdu
610041, P.R.China

Phone:86-28-85159085 Fax:86-28-85152372

website: www.aktinchem.com Email: info@aktinchem.com

CERTIFICATE OF ANALYSIS

Description	Oleuropein 98%min by HPLC		
CAS #	32619-42-4	MFW	540.5138g/mol
MF	C25H32O13	Manufacture Date	Nov.20,2018
Batch #	OL-181120	Quantity	13.8g

ITEMS	SPECIFICATION	RESULTS	METHODS
Appearance	White Powder	Conform	Visual
Identification	Positive	Positive	TLC&NMR
Oleuropein	≥98%	98.25%	HPLC
Loss on Drying	≤5%	1.03%	2H 105°C
Ashes	≤0.2%	Conform	2H 525°C

Storage: store in cool & dry place, keep away from strong light and heat

Stability: 24 months when properly stored

Packing: 1000mg/vial

Quality Control Manager: Ashily Zhou Quality Assurance Manager: Will Zhang