



Aktin Chemicals, Inc.

No.52, Yongfeng Road, Hi-tech Zone, Chengdu
610041, P.R.China

Phone:86-28-85159085 Fax:86-28-85152372

website: www.aktinchem.com Email: info@aktinchem.com

CERTIFICATE OF ANALYSIS

Description	Neochlorogenic acid		
CAS #	906-33-2	MFW	354.31g/mol
MF	C16H18O9	Manufacture Date	Dec.10,2018
Batch #	NCA-181210	Quantity	10.25g

ITEMS	SPECIFICATION	RESULTS	METHODS
Appearance	White Powder	Conform	Visual
Identification	Positive	Positive	TLC&NMR
Neochlorogenic acid	≥98%	99.52%	HPLC
Loss on Drying	≤5%	1.37%	2H 105°C
Ashes	≤0.2%	Conform	2H 525°C

Storage: store in cool & dry place, keep away from strong light and heat

Stability: 24 months when properly stored

Packing: 1000mg/vial

Quality Control Manager: Ashily Zhou Quality Assurance Manager: Will Zhang