



# Aktin Chemicals, Inc.

No.52, Yongfeng Road, Hi-tech Zone, Chengdu

610041, P.R.China

Phone:86-28-85159085 Fax:86-28-85152372

website: www.aktinchem.com Email: info@aktinchem.com

## CERTIFICATE OF ANALYSIS

<b>Description</b>	6-Gingerol 98%min by HPLC		
<b>CAS #</b>	23513-14-6	<b>MFW</b>	294.3859 g/mol
<b>MF</b>	C17H26O4	<b>Manufacture Date</b>	Dec.09,2016
<b>Batch #</b>	6G-161209	<b>Quantity</b>	18.5g

ITEMS	SPECIFICATION	RESULTS	METHODS
<b>Appearance</b>	Yellow Liquid	Conform	Visual
<b>Identification</b>	Positive	Positive	TLC&NMR
<b>6-Gingerol</b>	≥98%	98.29%	HPLC
<b>Loss on Drying</b>	≤5%	0.12%	2H 105°C
<b>Ashes</b>	≤0.2%	Conform	2H 525°C

**Storage:** store in cool & dry place, keep away from strong light and heat

**Stability:** 24 months when properly stored

**Packing:** 1000mg/vial

Quality Control Manager: Ashily Zhou Quality Assurance Manager: Will Zhang