



Aktin Chemicals, Inc.

No.52, Yongfeng Road, Hi-tech Zone, Chengdu
610041, P.R.China

Phone:86-28-85159085 Fax:86-28-85152372

website: www.aktinchem.com Email: info@aktinchem.com

CERTIFICATE OF ANALYSIS

Description	Procyanidin B2		
CAS #	29106-49-8	MFW	578.142426 g/mol
MF	C30H26O12	Manufacture Date	Aug.31,2017
Batch #	PR-170831	Quantity	2.27g

ITEMS	SPECIFICATION	RESULTS	METHODS
Appearance	Off-white powder	Conform	Visual
Identification	Positive	Positive	TLC&NMR
Procyanidin B2	≥98%	99.26%	HPLC
Loss on Drying	≤5%	1.14%	2H 105°C
Ashes	≤0.2%	Conform	2H 525°C

Storage: store in cool & dry place, keep away from strong light and heat

Stability: 24 months when properly stored

Packing: 1000mg/vial

Quality Control Manager: Ashily Zhou Quality Assurance Manager: Will Zhang